

## METHOD SUMMARY – QWI- FM0018

<b>Method Title</b>	Enterobacteriaceae - Petrifilm		
<b>Document number</b>	<b>QWI – FM0018</b>	<b>Date of Issued</b>	22 <sup>nd</sup> June 2017
Method External References	<p>3M Petrifilm™ Product Insert – Enterobacteriaceae Count</p> <p>Based on ISO 7402 1993 (Microbiology - General guidance for the enumeration of Enterobacteriaceae without resuscitation -- MPN technique and colony-count technique)</p> <p>Based on AS 5013.11.1-2004 (Food Microbiology - Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions of microbiological examination – General rules for the preparation of the initial suspension and decimal dilution)</p> <p>AFNOR Validated Method 3M 01/06-09/97</p> <p>AOAC Validated Method AOAC 2003.01-2006. Enumeration of Enterobacteriaceae in Selected Foods - Petrifilm Enterobacteriaceae Count Plate Method</p>		
Matrix	As listed on NATA Scope.		
ALS Department	<input type="checkbox"/> Pharmaceutical Chemistry <input type="checkbox"/> Water Microbiology <input checked="" type="checkbox"/> Food Microbiology <input type="checkbox"/> Pharmaceutical Microbiology <input type="checkbox"/> Food Chemistry		
Accreditation Status	<input checked="" type="checkbox"/> NATA (RYD) <input checked="" type="checkbox"/> NON-NATA (Other sites) <input type="checkbox"/> N/A		
Analysis technique	<input type="checkbox"/> HPLC <input type="checkbox"/> GC <input type="checkbox"/> Wet Chemistry <input type="checkbox"/> Physical <input type="checkbox"/> Gravimetric <input type="checkbox"/> Qualitative <input type="checkbox"/> Pour Plate <input type="checkbox"/> Spread Plate <input type="checkbox"/> MPN <input type="checkbox"/> Filtration <input checked="" type="checkbox"/> Petrifilm <input type="checkbox"/> EHS <input type="checkbox"/> ELISA <input type="checkbox"/> VIDAS UP <input type="checkbox"/> VIDAS <input type="checkbox"/> TEMPO		
Method Principle	<p>This method describes the procedure for the enumeration of <i>Enterobacteriaceae</i> by using a diagnostic rapid method designed for all foods and environmental samples (swabs). Petrifilm™ plates have not been tested with all possible food products, food processes testing protocols or with all possible strains of bacteria. This method is not suitable for Meat and Bone Meal samples, testing water, pharmaceutical or cosmetics.</p> <p>To check food samples for quality of production and contamination by potential pathogens. Specifically, <i>Enterobacteriaceae</i> are used to indicate the potential for enteric type organisms such as Salmonella to be present. It is an indicator of hygiene (post processing contamination) and a check for efficacy of processing (ie adequacy of cooking/heat treatment).</p> <p>This method can be used where a count of greater than 10 cfu/mL/g is expected in the sample. Validation of test is required before the use of this system.</p> <p>This method uses bacterial culture plates of dry medium and cold H<sub>2</sub>O soluble gel. Diluted or undiluted test suspensions are added to the plate at a rate of 1.0 mL per plate. Pressure, when applied to the plastic spreader</p>		

	<p>placed on overlay film, spreads test portions over 20 cm<sup>2</sup> growth area. The gelling agent is allowed to solidify and plates are incubated at 37 °C ± 1 °C for 24 ± 2 hours and then counted.</p> <p>An indicator in the Petrifilm™ <i>Enterobacteriaceae</i> count plate colours all colonies red. The top film traps gas produced by some bacteria. Acid-producing bacteria are seen as red colonies surrounded by yellow zones.</p> <p>Bacteria producing gas and/or acid are considered to be presumptive <i>Enterobacteriaceae</i>.</p>
Reporting Unit	Determination of Enterobacteriaceae cfu/g or cfu/mL
LOR/LOQ	<10 cfu/g or <10 cfu/mL

Minimum amount of sample required for analysis	20 g or 20 mL	Turnaround time	24 hours
--	---------------	-----------------	----------

<b>Author:</b>	Document Controller	<b>Date:</b>	20 <sup>th</sup> September 2017
<b>Authorised By:</b>	National Quality Manager	<b>Date</b>	20 <sup>th</sup> September 2017

CONFIDENTIAL