



	staphylococcus spp. is indicated by black colonies on BP due to the reduction of potassium tellurite and an egg yolk reaction. Typical and/or atypical colonies are confirmed by a Coagulase reaction, and the results given as the “presence” or “absence” of Coagulase positive staphylococcus spp.in x g or x mL of product.
Reporting Unit	Coagulase Positive Staphylococci Detected or Not Detected in 0.1 g or 1 g, 0.1 mL or 1 mL, or in amount tested as per Client request.
LOR/LOQ	Presence Absence (<1)

Minimum amount of sample required for analysis	10 g	Turnaround time	96 hours
--	------	-----------------	----------

<b>Author:</b>	Document Controller	<b>Date:</b>	21 <sup>st</sup> November 2017
<b>Authorised By:</b>	National Quality Manager	<b>Date</b>	21 <sup>st</sup> November 2017

CONFIDENTIAL